

**ELVE  
FLOW** an **ELVESYS** brand

# MUX RECIRCULATION VALVE 6/2

6-ports/2-positions

DOCUMENT REF: DTSMR6 201120

**DATASHEET**





## Introduction

The MUX Recirculation valve is a 6-ports / 2-positions microfluidic valve that is ideal for precise sample injection and fluid recirculation. When paired with the OB1 MKIII+, it can perform stable and unidirectional fluid recirculation or precise injection of a controlled volume of sample.

This versatile valve makes long-term experiments easier to automate and run.

## Main Features & Benefits

- Switch between 2 flow configurations in less than 180ms, thanks to the use of a fast & precise valve motor
- Work in total confidence thanks to the robust and biocompatible materials used (PTFE, PCTFE)
- No dead volume<sup>1</sup> and low carryover volume<sup>2</sup> (0.5µL), for the best solution usage efficiency.
- Easy 1/4-28 connection enables for a leak-free operation and quick connection using standard fittings.




<sup>1</sup> Volume that is stuck in the system (dead end), which is not clearly swept and relies on diffusion to clear out.


<sup>2</sup> Volume of liquid that will be mixed with the next liquid. It is not stuck, but will be swept next time a liquid passes.

# Product Specifications

## Description



**Figure 1 Front view of the MUX Recirculation:** The outer ports positions are marked from 1 to 6. Additional marking allows for the quick and easy identification of the 2 switching positions (A and B).



**Figure 2. Side view of the MUX Recirculation.** The USB and power supply connection ports can be found on the side of the device, along with the power button.

## Technical Specifications

<b>Performances</b>	<b>Port to port switching time (ms)</b>	180 ms
	<b>Max recommended pressure</b>	7 bar
	<b>Internal diameter</b>	0.5 mm
	<b>Internal volume <sup>3</sup></b>	2.5µL
	<b>Carryover volume <sup>4</sup></b>	0.5 µL
	<b>Dead volume <sup>5</sup></b>	None
	<b>Wetted materials</b>	PCTFE, PTFE
	<b>Number of ports</b>	6
	<b>Number of positions</b>	2
<b>Electrical</b>	<b>Input voltage range</b>	100V to 240V
	<b>AC supply frequency</b>	50 Hz to 60Hz
	<b>Power supply voltage</b>	18-24V DC
	<b>Max current consumption</b>	2A peak
	<b>Power consumption (max)</b>	36W
	<b>Connection type</b>	USB
<b>Mechanical specifications</b>	<b>Fluidic connectors</b>	Standard 1/4-28 UNF, flat-bottom
	<b>Operating temperature</b>	5-40°C
	<b>Operating Humidity</b>	20-70% non condensing
	<b>Dimensions (without connectors)</b>	133x156x133 mm
	<b>Weight</b>	1.12 kg
<b>Software</b>	<b>Computer specifications</b>	USB 2.0 port, Intel Pentium II 500 MHz, 1 Go Hard Disk space, 2 Go RAM Windows XP and newer, 2/64 bit. - LabVIEW® 2011 is required when using LabVIEW® libraries. - ESI installation is required to use SDK.
	<b>Provided software</b>	Elveflow Software Interface
	<b>Software development kit</b>	C++, Python, MATLAB® and LabVIEW® libraries

<sup>3</sup> Volume inside the system from entrance to exit.

<sup>4</sup> Volume of liquid that will be mixed with the next liquid. It is not stuck, but will be swept next time a liquid passes.

<sup>5</sup> Volume that is stuck in the system (dead end), which is not clearly swept and relies on diffusion to clear out.

## Product pack contents

Before setting up your MUX Recirculation valve, please check the package contents to verify that you have received the items below:


- the instrument
- a USB cable
- a power supply unit;

**Note:** The user documentation can be accessed anytime online through the Elveflow Support Portal (<https://support.elveflow.com/support/solutions>).

In addition to the above items, the user should have the necessary fluidic accessories, tubings and fittings to connect the inlets/outlets to the rest of the setup.

## Applications example

### Fluid recirculation with an Elveflow MUX Recirculation valve



**Illustration and setup principle to perform fluid recirculation with an Elveflow MUX Recirculation valve.** The detailed protocol for this application can be found in the MUX Recirculation valve user guide.

Other possible applications include (not limited to):

- [Controlled microfluidic sample injection](#)
- Fluid recirculation
- Cell culture on chip
- Drug screening
- Sample injection
- Stem cells assays
- Organ on chip
- microfluidics SPR/TIR imaging
- Toxicity tests