

Gx00-1L Industrial Peristaltic Pump Operating Manual



www.longerpump.com

Contents

1	Genera	al Information	3		
	1.1 Pr	ecaution	3		
	1.2 Re	epair Notes	3		
	1.3 Co	ontacts Information	3		
2	Produc	t Introduction	4		
	2.1 Main	Features	4		
	2.2 Unpa	cking	4		
	2.3 Syste	m Structure and Dimensions	4		
	2.4 Applie	cable Pump Head and Tubing	5		
	2.5 Syste	m Installation	8		
	2.5.1 P	Pump Head Installation	8		
	2.5.2 N	lounting of G100-1L/G300-1L/G600-1L			
	2.5.3 G	Grounding of G100-1L/G300-1L/G600-1L	13		
	2.6 Techr	nical Specifications	14		
	2.6.1 V	Vork Mode Introduction	14		
	2.6.2 C	15			
3	System	n Operation			
	3.1 Ke	eypad and Display			
	3.1.1 K	eypad Functions			
	3.1.2 D	Display Description			
	3.2 Sc	creens and Parameter Setting			
	3.2.1	Diagram of the Screen Changing Logic			
	3.2.2	Run Screens			
	3.2.3	Parameter Setting by Keypad			
	3.3 Flo	ow Rate Calibration			
	3.3.1	Flow Rate Calibration in Keypad Control Mode			
	3.3.2 F				
4	Definiti	on of the External Control Interface			
5	Append	dix			
	Appendix	1: Parameter Default Setting			
	Appendix	2: Modbus RTU Protocol			
	Appendix	3: Longer Pump OEM Communication Protocol	41		
	Appendix 4: Alarm Display				

1 General Information

1.1 Precaution

- The housing of the pump drive is equipped with a grounding terminal, in order to ensure the safety of the operator and improve the electromagnetic compatibility of the equipment, please reliably ground this product.
- The pump drive has high precision and high IP rating. Please do not disassemble, refit or maintain the product without permission, so as to avoid any damage or deterioration to its performance.
- All ports at the end of the drive shall be protected against damp and water. Besides, the end caps or plugs shall be screwed up in time.
- If tubing leaks or bursts, fluid may spray from the tubing and pump head. Take reasonable practicable measures to ensure the operator's safety.
- Make sure fluid in the tubing has been drained out, no pressure in the pipeline and disconnect pump from mains power, while removing or replacing the tubing.
- Disconnect the pump from the mains power before connecting the control signal wire.
- Do not touch the rotor while pump is running.
- Release the compression block when pump stop running for a long time to avoid tubing deformation caused by squeezing.
- Keep the rotor clean and dry to avoid tubing excessive wear and premature failure of pump head or drive.
- Please do not add the lubricating oil to the rotor by yourself, any improper operation could corrode the pump head housing or dislocate the tubing.
- Please connect the power cord, control cable, communication cable in correct way, and do not damage the plug.
- If the pump head can't resist organic solvent and corrosive liquid, please clean the liquid left on the surface of the pump head immediately

1.2 Repair Notes

Please contact Longer Pump or its distributor, and provide the product serial number before returning the product. Products which has been contaminated with, or exposed to, toxic chemicals or any other substance hazardous to health must be decontaminated before returning to Longer Pump or its distributor. Ship the product in its original packaging or better, to against possible damage or loss during the transport.

1.3 Contacts Information

Baoding Longer Precision Pump Co., Ltd 3rd Floor, Building 6B, University Science Park Baoding National, High - Tech Industrial Development Zone, Baoding, Hebei, China 071051 Email: <u>longer@longerpump.com</u> Tel: 86-312-3110087 Fax: 86-312-3168553 www.longerpump.com

2 Product Introduction

The Gx00-1L series is an industrial peristaltic pump drive specially developed for severe industrial applications. It is small and simple in shape and features a compact structure. With good electromagnetic compatibility, the pump will not fail in the EMI environment and will not interfere with other electronic equipment nearby. It also has high IP rating (IP66), which is suitable for wet and dusty environment. In addition, it has aluminum alloy housing with good heat dissipation performance, and the anti-corrosion paint treatment improves the corrosion-resistant performance. The pump drives are compatible with many Longer Pump peristaltic pump heads and accepts many different tubings. G100-1L can provide a max flow rate of 750mL/min (single channel), G300-1L can provide a max flow rate of 1.75L/min (single channel), G600-1L can provide a max flow rate of 3.25L/min.

2.1 Main Features

The main functions of Gx00-1L series products are as follows:

- Support multiple work modes: continuous transfer, timed transfer, volume transfer, prime transfer at full speed
- Calibration function: calibrate the flow rate through input actual flow rate or pumped volume
- LCD can display a variety of parameters: pump speed, flow rate, running time and pumped volume can be displayed on different screens.
- Support three control modes: keypad control, RS485 communication control, and external signal control.

2.2 Unpacking

- (1) Take out the pump and accessories from the shipping carton.
- (2) Check the packing list and confirm that the attachments are intact and complete.
- (3) In case of any questions, please contact Longer Pump or local distributor.

2.3 System Structure and Dimensions

The Gx00-1L series peristaltic pump includes the following parts:

- ① Pump Drive
- ② Pump Head + Tubing





Outline Dimension

Fuse: F2A 250V

2.4 Applicable Pump Head and Tubing

Pump Drive Model	Pump Head	Applicable Silicone Tubing	Applicable Pharmed Tubing	Max Flow Rate Reference mL/min
	YZ1515x,	13#, 14#, 19#, 16#, 25#,	13#, 14#, 19#, 16#, 25#,	570
	YZ II 15	17#, 18#	17#, 18#	570
	FG15-13	13#, 14#, 19#, 16#, 25#, 17#, 18#	13#, 14#, 19#, 16#	645
G100-1L	DMD15-13-B	2*13#, 2*14#, 2*19#, 2*16#, 2*25#	2*13#, 2*14#, 2*19#, 2*16#	560
	YZ2515x	15#, 24#	Not recommended	400
	YZ II 25	15#, 24#, 35#, 36#	Not recommended	750
	FG25	15#, 24#	Not recommended	480
	BZ15-13-A	14#	14#	33

BZ15-13-B	16#	16#	120	
BZ15-13-C	25#	25#	225	
BZ15-13-D	17#	17#	400	
BZ25-13-B	24#	Not recommended	375	
DG15-24	16#, 25#, 17#	Not recommended	450 (single channel)	
DG15-28	13#, 14#, ID≤3.17mm, wall thickness 0.8-1mm	ID≤3.17mm Wall thickness: 0.8-1mm	72 (single channel, pump speed ≤100rpm)	
DG-1-A (C)	ID≤3.17mm Wall thickness: 0.8-1mm	ID≤3.17mm Wall thickness: 0.8-1mm		
DG-2-A (C)	ID≤3.17mm Wall thickness: 0.8-1mm	ID≤3.17mm Wall thickness: 0.8-1mm		
DG-4-A (C)	ID≤3.17mm Wall thickness:0.8-1mm	Not recommended	48 (single channel, pump speed ≤100rpm)	
DG-6-A (C)	ID≤3.17mm Wall thickness: 0.8-1mm	Not recommended		
DG-8-A (C)	ID≤3.17mm Wall thickness: 0.8-1mm			
DG-1-B (D)	ID≤3.17mm Wall thickness: 0.8-1mm	ID ≤3.17mm Wall thickness: 0.8-1mm		
DG-2-B (D)	ID ≤3.17mm Wall thickness: 0.8-1mm	ID≤3.17mm Wall thickness: 0.8-1mm	32 (single channel, pump speed ≤100rpm)	
DG-4-B (D)	ID≤3.17mm Wall thickness: 0.8-1mm	Not recommended		
dPOFLEX BPH01	13#, 14#, 19#, 16#, 25#		130	

Pump Drive Model	Pump Head	Applicable Silicone Tubing	Applicable Pharmed Tubing	Max Flow Rate Reference mL/min
	YZ1515x, YZ II 15	13#, 14#, 19#, 16#, 25#, 17#, 18#	13#, 14#, 19#, 16#, 25#, 17#, 18#	1280
	FG15-13	13#, 14#, 19#, 16#, 25#, 17#, 18#	13#, 14#, 19#, 16#	1400
G300-1L	DMD15-13-B	2*13#, 2*14#, 2*19#, 2*16#, 2*25#	2*13#, 2*14#, 2*19#, 2*16#	1200
	YZ2515x	15#, 24#	Not recommended	930
	YZ 25	15#, 24#, 35#, 36#	Not recommended	1750
	FG25-13	15#, 24#	Not recommended	1280
	BZ15-13-A	14#	14#	88

BZ15-13-B	16#	16#	265
BZ15-13-C	25#	25#	560
BZ15-13-D	17#	17#	930
BZ25-13-B	24#	Not recommended	930
DG15-24	16#, 25#, 17#	Not recommended	1050 (Single Channel)
BPH01	13#, 14#, 19#, 16#, 25#	13#, 14#, 19#, 16#, 25#	300

Pump Drive Model	Pump Head	Applicable Silicone Tubing	Applicable Pharmed Tubing	Max Flow Rate Reference mL/min
	YZ1515x, YZ II 15	13#, 14#, 19#, 16#, 25#, 17#, 18#	13#, 14#, 19#, 16#, 25#, 17#, 18#	2380
	FG15-13	13#, 14#, 19#, 16#, 25#, 17#, 18#	13#, 14#, 19#, 16#	2600
	DMD15-13-B	2*13#, 2*14#, 2*19#, 2*16#, 2*25#	2*13#, 2*14#, 2*19#, 2*16#	2500
	YZ2515x	15#, 24#	Not recommended	1730
	YZ II 25	15#, 24#, 35#, 36#	Not recommended	3250
	FG25-13	15#, 24#	Not recommended	2380
G600-1L	BZ15-13-A	14#	14#	160
	BZ15-13-B	16#	16#	490
	BZ15-13-C	25#	25#	1040
	BZ15-13-D	17#	17#	1730
	BZ25-13-B	24#	Not recommended	1730
	DG15-24	16#, 25#, 17#	Not recommended	1900
	BPH01 (intermittent when speed>400rpm)	13#, 14#, 19#, 16#, 25#	13#, 14#, 19#, 16#, 25#	530 (@600rpm)

2.5 System Installation

2.5.1 Pump Head Installation

(1) Installation of Pump Head YZ1515x/YZ2515x/YZII15/YZII25

Insert the tang of the pump head shaft into the slot of the motor shaft, and align the two mounting holes on the pump head with the mounting holes on the pump drive. Then insert the two mounting screws into the mounting holes and tighten the screws.



(2) Installation of Pump Head FG15-13/FG25-13

Mount the pump head:

Mount the mounting plate on the drive unit by tightening three mounting screws M4X10. Insert the tang of the pump head shaft into the slot of the motor shaft, turn the pump head 45 degrees relative to the vertical direction, engage the bayonet on the mounting plate with the bayonet slot on the back of the pump head, turn the pump head clockwise until it locks into an upright position.

Remove the pump head:

Push the locking lever back and turn the pump head anti-clockwise about 45 degrees. Then take off the pump head from the mounting plate.



(3) Installation of Pump Head DMD15-13 and Tubing

Mount the pump head:

Mount the mounting plate on the drive unit through three cross recess head screws M4X10. Release the two levers to open the pump head and remove the compression block. Insert the tang of the pump head shaft into the slot of the motor shaft. Press the pump head firmly against the mounting plate. Turn the pump head to align the mounting holes on the pump head with the mounting holes on the mounting plate. Insert the two mounting screws (hexagon socket head cap screw M3X8) into the mounting holes and tighten the screws.



Load the tubing:

Release the levers to remove the compression block. Insert the tubing fitting assembly into the compression block.

Note: When using 25# tubing, move the partition blocks to the outermost positioning hole.



Put the compression block with tubing fitting assembly back to the pump head, and lock the levers.



(4) Installation of Pump Head BZ15/BZ25

Insert the tang of the pump head shaft into the slot of the motor shaft, and align the four mounting holes on the pump head with the mounting holes on the pump drive. Then insert the four mounting screws into the mounting holes and tighten the screws.



(5) Installation of Pump Head DG15-24/DG15-28

Mount the mounting plate on the drive unit through four screws M4X10. Remove the tubing and the pump head cartridges. Insert the tang of the pump head shaft into the slot of the motor shaft. Press the pump head firmly against the pump drive. Turn the pump head to align the mounting holes on the pump head with the mounting holes on the pump drive. Insert the screws M4X16 and tighten the screws.



(6) Installation of Pump Head DG Series

Remove the tubing and the pump head cartridges. Insert the tang of the pump head shaft into the slot of the motor shaft. Press the pump head firmly against the pump drive. Turn the pump head to align the mounting holes on the pump head with the mounting holes on the pump drive. Insert the 2 hexagon socket head cap screws M4X8 and tighten the screws.





Note: The adapter bracket is an optional accessory.



G100-1L, G300-1L and G600-1L are all equipped with grounding rod. Reliable grounding of the pump will prevent the misoperation and damage from electromagnetic interference.

Grounding steps:



Unscrew the grounding knob on the rear of the drive. Do not loosen the part fixed on the housing, as shown in above diagram. Insert the grounding cable through the hole on the grounding rod

Tighten the grounding knob to fix the grounding cable on the grounding rod.

Or directly solder the grounding cable on the grounding rod.



2.6 Technical Specifications

Pump model	G100-1L	G300-1L	G600-1L		
Speed range	0, 0.01-150 (rpm), CW/CCW	0, 0.01-350 (rpm), CW/CCW	0, 0.01-650 (rpm), CW/CCW		
		0.01rpm (< 10rpm)			
	0.	1rpm (10rpm≤speed < 100rpm)			
Speed resolution	1rpi	m (100rpm≤speed≤Max speed)			
Flow rate range	0.15uL/min-750mL/min	0.15uL/min-1750mL/min	0.15uL/min-3250mL/min		
Display	LCD for speed, flow rate, run	ning time, pumped volume, and ot	her setting parameters		
Work mode	Continuous	transfer, timed transfer, volume tra	insfer		
Calibration	Improve the fl	ow rate and dispensing volume ac	curacy		
	Keypad control (start/stop can	be controlled by a footswitch), ext	ernal signal control and		
Control mode		communication control			
	Start/stop control, direction control: log	gic level signal or switch signal (dry	 contact), momentary or 		
	maintained trigger mode can be config	jured			
	Speed control: 0-5V, 0-10V, 4-20mA,	0-10kHz signal to control the pump	o speed		
	External control parameters can be configured through keypad (signal type, trigger mode, max speed),				
External control and communication commands (signal type, trigger mode and speed/signal range).					
	RS485 interface, Modbus RTU and Longer Pump OEM protocol, communication control paramet				
Communication control	can be configured through keypad and communication commands				
Keypad can be locked to prevent misoperation, the delay time (5s to 1min) before lock			o 1min) before lock can be		
Keypad lock	configured				
	Running parameters and system parameters can be saved automatically (it takes 1 second to save the				
Parameter memory	parameters)				
Pump status when					
powered up	The pump status when powered up can be set to stop or the same status before power-off				
Prime	Fast filling or emptying at full speed				
	CE, UKCA				
	EMC: EN IEC61000-6-2:2019 EN IEC61000-6-4:2019				
Certificate	LVD: EN 61010-1:2010/A1:2019 EN 61010-1:2010/A1:2019				
Dimension (LxWxH)	240*142*152(mm)				
Power supply	upply AC100-240V 50Hz/60Hz				
Power consumption	25W	40W	60W		
Working condition	Environment temperature: 0° C ~ 40° C, relative humidity: $10\% \sim 90\%$				
IP rating		IP66			
Weight	3.64kg	3.75kg	3.75kg		

2.6.1 Work Mode Introduction

Gx00-1L series pump can work in different work modes: continuous transfer mode, timed transfer mode and volume transfer mode.

Continuous transfer mode: the pump will run continuously and will stop when receiving a stop signal from the keypad, external signal or communication command.

Timed transfer mode: set the target time before starting the pump, the pump will begin timing when starts running, and the pump will stop automatically when the target time is reached.

Volume transfer mode: set the target volume before starting the pump, the pump will begin timing when starts running, and the pump will stop automatically when the target volume is reached. (volume = flow rate * time). Before working in the volume transfer mode, the flow rate should be calibrated so that the actual pumped volume can match the target volume.

2.6.2 Control Mode Introduction

There are three control modes: keypad control, external signal control, communication control. Refer to Chapter 3.3.3 for control mode setting details.

(1) **Keypad control mode (INT)**: set all the parameters and control the pump running through the keypad.

In the keypad control mode, the start/stop of the pump can also be controlled by a logic level signal, a switch signal (dry contact signal), or a footswitch. Refer to Chapter 3.3.3 for the start/stop trigger mode setting details.

Note:

- 1. The keypad control mode can only be set through the keypad.
- **2.** In the keypad control mode, the Start/stop key, a footswitch, or other external start/stop signal are all valid.
- (2) **External signal control mode (EXT)**: control the start/stop status, running direction and pump speed by external signals.

Before controlling the pump with external signals, the parameters of speed signal type, max speed corresponding to max speed signal, start/stop trigger mode, and run direction trigger mode should be set through the keypad or communication command. Refer to Chapter 3.3.3 for parameter setting by keypad. Refer to Appendix B for parameter setting by command. Refer to Chapter 4 for the signal requirement and wiring instructions.

Pump speed: can be controlled by 0-5V, 0-10V, 4-20mA, 0-10kHz signal, and max speed corresponding to the max signal can be set. The input speed signal should be the same as the signal type setting. when the pump is working in the external signal control mode, tap Enter Key to check the current speed signal setting.

Start/stop: can be controlled by a logic level signal (5V-24V) or a switch signal (dry contact signal). The trigger mode can be set to pulse trigger (maintained) or level trigger (momentary).

Running direction: can be controlled by a logic level signal (5V-24V) or a switch signal (dry contact signal). The trigger mode can be set to pulse trigger (maintained) or level trigger (momentary).

Note:

- 1. The external control mode only can be set through the keypad.
- 2. In the external signal control mode, the start/stop and running direction can not be controlled through the keypad. But press and hold the Start/stop key, the pump will stop.
- 3. In the external signal control mode, the timed transfer work mode and volume transfer work mode are invalid.
- (3) **Communication control mode (COM)**: set all parameters and control the pump running through RS485 interface and communication commands. (Both Modbus RTU and Longer Pump OEM protocol are supported.)

Before controlling the pump with communication commands, the parameters of pump address, baud rate and parity should be set through the keypad or communication command. Refer to Chapter 3.3.3

for parameter setting by keypad. Refer to Appendix B for parameter setting by command. Refer to Chapter 4 for the wiring instructions.

Note:

- 1. The communication control mode can be set through the keypad or communication command.
- 2. In the communication control mode, the start/stop and running direction can not be controlled through the keypad. But press and hold the Start/stop key, the pump will stop.
- 3. Communication control has the highest priority. When the pump is working in keypad control or external signal control mode, send the communication control enable command, the pump will stop running and change to communication control mode.
- 4. In the communication control mode, the flow rate can be calibrated by send the actual pumped volume or actual flow rate value. But the calibration screen will not be displayed on the screen.

3 System Operation

3.1 Keypad and Display



3.1.1 Keypad Functions

Кеу	Functions	Remark	
		The key is invalid in the	
		external signal control	
\frown	In the keypad control mode, when the pump is on	mode, communication	
26	standby with a running screen or is working in the	control mode, timed	
Bunning direction key	continuous transfer mode, tap this key to switch the	transfer mode, volumed	
Running direction key	running direction.	transfer mode, calibration	
		process, and parameter	
		setting screen.	
	1. Change the decimal point position: when the pump		
	is on standby with a running screen, tap the unit switch		
	key to change the decimal point position and units of		
Unit switch key	the pump speed, flow rate, and dispensing volume.		
	2. Change the calibration variable: when the pump		

	isplays the calibra			
	o switch the varia			
	ctual flow rate to			
	ctual pumped volu	ime to calibrate).		
	. Set the value:	when the pump displays a running		
	screen, tap tl	ne key "+" or "-" to increase or		
	decrease the	owest digit of the displayed value.	When the pump is running,	
	Press and ho	d the keys to quickly increase or	the "+" and "-" keys are only	
	decrease the c	lisplayed value.	valid to change the pump	
Increase and decrease keys	. Change the p	arameter options: when the pump	speed and flow rate.	
	displays the pa	rameter setting screen, tap the key		
	"+" or "-" to cha	ange the parameter options.		
	. In the keypad	control mode: when the pump		
	displays a run	ning screen, tap the key to start or	In the keypad control mode,	
	stop the pump		when the pump is	
	. In the externa	control mode and communication	connected to a footswitch or	
Start/ stop key	control mode,	when the pump is running, press	other external start/stop	
	and hold the s	tart/stop key to stop the pump. The	signal, the Start/stop key is	
	control mode v	vill not be changed.	still valid.	
	. Full speed: In	the keypad control mode, when the		
	pump dislays a	running screen, press and hold the	1. Full speed function is	
	Full speed key	to run the pump at full speed, and	invalid in timed	
	release the ke	y to return to the previous status.	transfer mode,	
	. Confirm the pa	arameter: when the pump displays	volumed transfer	
	the parameter	setting screen, tap the Enter key to	mode and calibration	
	confirm the se	tting value and proceed to the next	process.	
Full speed/Enter key	parameter sett	ing screen.	2. The setting	
	. Pump calibrat	ion: when the pump displays the	parameters in the	
	calibration scr	een, tap the Enter key to calibrate	parameter screen will	
	the pump acco	ording to the displayed flow rate or	be invalid without	
	volume.		tapping the Enter key.	
	. Change the	unning screen: when the pump		
	displays the ru	inning screen, tap Function key to		
	change the r	unning screen. When the pump		
	displays the	parameter setting screen, tap		
	Function key t	o cancel current parameter setting		
	and return to p	revious parameter setting screen.		
	. Enter/exit the p	parameter setting screen: press and		
	hold the Fur	ction key to enter or exit the		
	parameter set	ting screen. The first parameter		
Function Key	setting screen	is the control mode setting screen.		
	Exit action car	be done when pump displays any		
	parameter sett	ing screen. And this function is valid		
	in any control ı	node.		
	. Unlock the ke	ypad: when the keypad is locked,		
	press and hol	d the Function key to unlock the		
	keypad.			

3.1.2 Display Description



lcon	Name	Description	Remark
	Run	This icon is displayed on the running screen when the	
		pump is running	
	Ston	This icon is displayed on the running screen when the	
	5100	pump is stopped	
5	Counter-	This icon is displayed on the running screen when the	
	clockwise	pump is running counter-clockwise.	
	Cleakwise	This icon is displayed on the running screen when the	
51	CIOCKWISE	pump is running clockwise.	
		This ison annears on the maning senser when the	The keypad lock function can be enabled in
\bigcirc	Keypad lock	This icon appears on the running screen when the	the parameter setting screen, and the delay
		кеуюск із юскей	time (5s to 1min) before lock can be set.
	Keypad control	This icon is displayed on the running screen when the	
I IN I	mode	current control mode is keypad control.	
EVT	External control	This icon is displayed on the running screen when the	
EXI	mode	current control mode is external control.	
0014	Communication	This icon is displayed on the running screen when the	
GUI	control mode	current control mode is communication control.	
0.41		This icon on the running screen indicates that the flow	
UAL	Calibration	rate can be calibrated on the current screen.	
M-M-M	Parameter and	This area displays pump system parameters, running	
1:1:1	status display	parameters and pump status.	
		, IIII and indicate the current value is	
		ul /m	
uL/m	Unit for volume	volume.	
mL/m	or flow rate	···· /m ···· ·· ·· /m ··· /m	
		mL/min.	
		/m	
		Indicate the current value is flow rate.	
rnm	Unit for pump	This icon on the running screen indicates that the	
r pm	speed	current value is pump speed.	

3.2 Screens and Parameter Setting

3.2.1 Diagram of the Screen Changing Logic



3.2.2 Run Screens

Gx00-1L series pump can work in different work modes: continuous transfer mode, timed transfer mode and volume transfer mode. When the pump is running, the pump speed, flow rate, time and volume can be

displayed on the screen. 5 run screens can be displayed on the LCD to indicate the pump running status. Tap

the Function key 🤹 to switch the screens. The pump can start running when displaying any run screen.



Run Screen	Display Example	Description	
Pump Speed Screen	The information displayed: Keypad control mode, pump speed 600rpm, stop status, clockwise direction	On the Pump Speed Screen, the displayed indicating the three-digit number on the screed pump speed. The screen also shows the start/stop state direction and the control mode. When the pump is powered on for the first time display the pump speed screen by default. The pump speed can be adjusted when the pu- on standby. The pump speed and flow rate correspond to value will be automatically changed as the oth After calibration is completed, the pump spee match the set flow rate. If the calibrated pump the max speed, the pump speed will adjust to and the flow rate will change to the max flow speed. On the Flow Rate Screen, the displayed unit is	ed unit is rpm , in represents the tus, the running e, the screen will ump is running or each other. One her is adjusted. eed will adjust to o speed exceeds o the max speed, v rate at the max
Flow Rate Screen	The information displayed: Keypad control mode, flow rate 600mL/min, stop status, clockwise direction	 L/m, indicating the three-digit number represents the flow rate. The screen also shows the start/stop state direction and the control mode. The flow rate can be adjusted when the pump standby. The pump speed and flow rate correspond to value will be automatically changed as the oth The flow rate should be calibrated when runn the first time, after replacing the pump head, of the tubing. So that the displayed flow rate will flow rate. 	on the screen tus, the running o is running or on each other. One her is adjusted. ing the pump for or after replacing match the actual
Time Screen	The information displayed: Keypad control mode, set the target time to 1 minute, stop status, clockwise direction	On the Time Screen, the displayed three represents time. The screen also shows the start/stop state direction and the control mode. When the pump is on standby, the time care seconds to 999 minutes. When the time is set will run in continuous transfer mode. When t other values, the pump will run in timed transfer	ee-digit number tus, the running n be set from 0 : as 0s, the pump the time is set to fer mode.

		۶	When the pump starts running, it will begin timing. The screen
			will show the elapsed running time. When the preset target
			set to 0s, and the number running time exceeds the max display
			value of 999min the displayed time will start again from 0s
			1:00 means 1 minute, 10:1 means 10 minutes and 10 seconds
		-	100 means 100 minutes
		Δ	On the Volume screen, the displayed unit is up or the
		-	indicating the three digit number on the screen represents the
			volume
		2	The screen also shows the start/stop status the running
		Í	direction and the control mode
		Δ	When the nump is on standby the volume can be set from 0 to
		Í	Q001 When the volume is set to 0 the nump will run in
			continuous transfer mode. When the volume is set to other
			values the nump will run in volume transfer mode
	INT	2	When the nump starts running, it will begin calculating the
Volume Screen	The information displayed:	ĺ	numbed volume based on the flow rate displayed on the Flow
	Keypad control mode set the target		Rate Screen and the elansed running time. And the screen will
	volume to 600ml stop status		show the number volume. When the preset target volume is
	clockwise direction		reached the nump will automatically stop. If the volume is set
			to Oul, and the numbed volume exceeds the max display value
			of 9991 the displayed volume will start again from Oul
			The flow rate should be calibrated when running the pump for
			the first time, after replacing the pump head, or after replacing
			the tubing. So that the displayed volume will match the actual
			pumped volume.
	CAL CAL	>	Indicates the current screen is the Calibration Screen.
	The information displayed:		The Calibration Screen is used to calibrate the flow rate. After
	Keypad control mode, calibrate the		calibration, the displayed flow rate on the Flow Rate Screen
	flow rate with the pumped volume of		and the pumped volume on the Volume Screen will match the
	600mL, stop status, clockwise		actual values.
Calibration	direction	۶	The flow rate can be calibrated by using the actual pumped
Screen			volume or actual flow rate. Press and hold the Unit Switch Key
			to switch the unit between the volume unit and the flow
			rate unit. Refer to Chapter 3.4 for the calibration process
			details.
	I he information displayed:	>	I he flow rate should be calibrated when running the pump for
	Keypad control mode, calibrate the		the first time, after replacing the pump head, or after replacing
	tiow rate with the actual flow rate of		the tubing.
	600mL/m, stop status, clockwise		
	direction		

Note:

- 1. When the pump is powered up, the LCD will display one of the run screens: Pump Speed Screen, Time Screen, or Volume Screen, depending on the pump work mode before power off. If the pump was in the continuous work mode (target time and target volume were set to 0), the pump will display the Pump Speed Screen. If the pump was in the timed transfer mode (target time was set to a value other than 0, and the volume was calculated based on the time and flow rate), the pump will display the Time Screen. If the pump was in the volume transfer mode (target volume was set to a value other than 0, and the time was calculated based on the time and flow rate), the pump will display the Time Screen.
- 2. When the pump speed and flow rate are 0, the pump can not start.
- 3. When the target volume or target time is set too small, the pump can not start
- 4. When pump is working in the external signal control mode or communication control mode, the LCD can display Pump Speed Screen, Flow Rate Screen, Time Screen and Volume Screen.

3.2.3 Parameter Setting by Keypad

In any control mode and on any run screen, when the pump is on standby, press and hold the Function Key

to enter the Parameter Setting Screen. Refer to the below table for the parameters that can be set through the keypad and the options.

When the pump displays any parameter setting screen, tap the Function Key

to switch the setting screen for different parameters. Tap Increase and Decrease Keys

among the parameter options. After the parameter setting, tap the Enter Key to confirm and save the

setting, and show the next parameter setting screen. When all the parameters have been set up, press and



hold the Function Key 🔀 to return to the previous run screen, then the pump can be controlled to start.

ltem	Parameters	Options
1	Control Mode	INT: for keypad control (start/stop can be controlled by the keypad, a footswitch and an
		external signal)
		EXT: for external signal control
		COM: for communication control
2	Pump status when	
	powered up	when powered up, the pump will return to its previous status before power-off
	(only valid for INT mode	
	and continuous transfer	LI: when powered up, the pump will stop
	mode)	
3	Keypad lock	OFF: keypad will not lock
		0:05-1:00 (5s-1min): keypad will lock when the keypad is not operated within the set time
4	External control speed	0-5V, 0-10V, 4-20mA, 0-10kHz
	signal	
	(only valid for the EXT	

	mode)	
5	Max speed for the	Setting range for G100-1L: (min speed+1rpm) to 150rpm
	external control signal	Setting range for G300-1L: (min speed+1rpm) to 350rpm
	(only valid for the EXT	Setting range for G600-1L: (min speed+1rpm) to 650rpm
	mode)	The default value of min speed is 0rpm. The min speed can be set by communication
		command.
6	Start/Stop control trigger	Pulse trigger (maintained)
	mode	Level trigger (momentary)
7	Running direction control	Pulse trigger (maintained)
	trigger mode (only valid	Level trigger (momentary)
	for the EXT)	
8	Pump address for the	1-30
	communication control	
9	Baud rate for the	1200bps, 2400bps, 4800bps, 9600bps, 19200bps, 38400bps, 115200bps
	communication control	
10	Parity for the	None
	communication control	EVEN
		ODD



3.2.3.1 Control Mode Setting

There are three control modes: keypad control mode, external signal control mode and communication control mode.

In the keypad control mode, the start/stop can be controlled by the keypad, a footswitch and an external signal.

Setting steps:	
----------------	--

1. When the pump is on standby and displaying any run screen, press and hold the Function Key enter the Control Mode setting Screen.

External Signal Control

- 2. Tap Increase and Decrease Keys to switch among the control modes. INT means keypad control, EXT means external signal control, COM means communication control.
- to confirm and save the control mode setting, and display the next parameter 3. Tap the Enter Key setting screen (Pump Status When Powered Up).
- 4. If no other parameters need to be set, press and hold the Function Key to return to the previous run screen. The pump only can start running when displaying run screen.

3.2.3.2 Pump Status When Powered Up

When the control mode is the keypad control and work mode is continuous transfer, the pump status when powered up can be set to stop or the status before power-off.



When powered up, pump will stop

Setting steps:

1. When the pump is on standby and displaying any run screen, press and hold the Function Key

enter the Control Mode setting screen. Then tap Enter Kev to enter the Pump Status When Powered Up setting screen.

- 2. Tap Increase and Decrease Keys
- to switch between the two statuses.



up, pump will return to its previous status before power-off.





its previous status before power-off

Communication Control







when powered



- Tap the Enter Key to confirm and save the setting, and display the next parameter setting screen (Keypad Lock).
- 4. If no other parameters need to be set, press and hold the Function Key



to return to the previous

run screen. The pump only can start running when displaying run screen.

Note: the setting for pump status when powered up is only valid for the keypad control mode and the continuous work mode. When pump is working at timed mode or volume mode, the pump will always stop when powered up.

3.2.3.3 Keypad lock

When the pump displays the run screen, the keypad can lock automatically to prevent misoperation. The time before lock can be set to 5s - 1min.



The keypad will never lock



The keypad will lock when it is not operated within 1min

Setting steps:

1. When the pump is on standby and displaying any run screen, press and hold the Function Key \square

enter the Control Mode setting screen. Then tap Enter Key wice to enter the Keypad Lock setting screen.

- 2. Tap Increase and Decrease Keys to set the keypad lock delay time.
- 3. Tap the Enter Key to confirm and save the setting, and display the next parameter setting screen (Speed Signal).
- 4. If no other parameters need to be set, press and hold the Function Key 🐨 to return to the previous run screen. The pump only can start running when displaying run screen.

Note:

1. When the keypack is locked, press and hold the Function Key

^J to unlock the keypad.

2. Keypad lock function is valid for keypad control mode, external signal control mode and communication control mode.

3.2.3.4 External Control Speed Signal

When using an external signal to control the pump, the pump speed signal can be 4-20mA, 0-5V, 0-10V, 0-10kHz. Before controlling the pump, the speed signal should be set up on the pump.



Setting steps:

1. When the pump is on standby and displaying any run screen, press and hold the Function Key

enter the Control Mode setting screen. Then tap Enter Key three times to enter the Speed Signal setting screen.

- 2. Tap Increase and Decrease Keys to set the speed signal.
- 3. Tap the Enter Key to confirm and save the setting, and display the next parameter setting screen (Max Speed for External Control).
- 4. If no other parameters need to be set, press and hold the Function Key to return to the previous run screen. The pump only can start running when displaying run screen.

Note: when the pump is in the external signal control mode, and the pump is displaying the run screen, tap

Enter Key to check the speed signal setting.

3.2.3.5 Max Speed for the External Control

In the external signal control mode, the max speed corresponding to the max speed signal can be set. Setting range for G100-1L: (min speed+1rpm) to 150rpm Setting range for G300-1L: (min speed+1rpm) to 350rpm Setting range for G600-1L: (min speed+1rpm) to 650rpm The default value of min speed is 0rpm. The min speed can be set through communication command.



The max speed is set to 100rpm

Setting steps:

1. When the pump is on standby and displaying any run screen, press and hold the Function Key



- 2. Tap Increase and Decrease Keys to set the max speed corresponding to the max speed signal.
- to confirm and save the setting, and display the next parameter setting screen 3. Tap the Enter Kev (Start/stop Trigger Mode).
- to return to the previous 4. If no other parameters need to be set, press and hold the Function Key

run screen. The pump only can start running when displaying run screen.

3.2.3.6 Start/stop Control Trigger Mode

When using an external logic level signal or a switch signal (dry contact signal or footswitch) to start or stop the pump, the trigger mode can be set to pulse trigger (maintained) or level trigger (momentary). Through the keypad, the pulse trigger mode can be set to failing edge trigger, the level trigger mode can be set to high level signal trigger. Through the communication commands, the pulse trigger mode can be set to failing edge trigger or rising edge trigger, the level trigger can be set to high level signal trigger or low level signal trigger.



Failing Edge Trigger

Rising edge Trigger





Low Level Signal Trigger

Failing Edge Trigger: When the control signal is a logic level signal, a change from high (5V to 24V) to low level is considered a failing edge signal. When the control signal is a switch signal (dry contact), a change from open to closed contact is considered a failing edge signal. Upon receiving a failing edge signal, the pump will toggle its status between start and stop.

Rising Edge Trigger: When the control signal is a logic level signal, a change from low to high level (5V to 24V) is considered a rising edge signal. When the control signal is a switch signal (dry contact), a change from closed to open contact is considered a rising edge signal. Upon receiving a rising edge signal, the pump will toggle its status between start and stop.

High Level Signal Trigger: When the control signal is a logic level signal, a voltage of 5V to 24V is considered high level. When the control signal is a switch signal (dry contact), an open contact is considered high level, and a closed contact is considered low level. When the signal changes from low level to high level, and remains high level, the pump will start.

Low Level Signal Trigger: When the control signal is a logic level signal, a voltage of 5V to 24V is considered high level. When the control signal is a switch signal (dry contact), an open contact is considered high level, and a closed contact is considered low level. When the signal changes from high level to low level, and remains low level, the pump will start.

Setting steps:

1. When the pump is on standby and displaying any run screen, press and hold the Function Key

enter the Control Mode setting screen. Then tap Enter Kev five times to enter the Start/Stop Trigger Mode setting screen.

- 2. Tap Increase and Decrease Kevs to set the triager mode.
- to confirm and save the setting, and display the next parameter setting screen 3. Tap the Enter Kev (Running Direction Trigger Mode).
- 4. If no other parameters need to be set, press and hold the Function Key to return to the previous

run screen. The pump only can start running when displaying run screen.

3.2.3.7 Running Direction Control Trigger Mode

When using an external logic level signal or a switch signal (dry contact signal or footswitch) to control the running direction of the pump in the external signal control mode, the trigger mode can be set to pulse trigger (maintained) or level trigger (momentary). Through the keypad, the pulse trigger mode can be set to failing edge trigger, the level trigger mode can be set to high level trigger. Through the communication commands, the pulse trigger mode can be set to failing edge trigger or rising edge trigger, the level trigger mode can be set to high level signal trigger or low level signal trigger.

EXT	EXT	EXT	EXT

Failing Edge Trigger

Rising edge Trigger

High Level for CW



Failing Edge Trigger: When the control signal is a logic level signal, a change from high (5V to 24V) to low level is considered a failing edge signal. When the control signal is a switch signal (dry contact), a change from open to closed contact is considered a failing edge signal. Upon receiving a failing edge signal, the pump will toggle its running direction between clockwise and counterclockwise.

Rising edge Trigger: When the control signal is a logic level signal, a change from low to high level (5V to 24V) is considered a rising edge signal. When the control signal is a switch signal (dry contact), a change from closed to open contact is considered a rising edge signal. Upon receiving a rising edge signal, the pump will toggle its running direction between clockwise and counterclockwise.

High Level Signal for CW: When the control signal is a logic level signal, a voltage of 5V to 24V is considered high level. When the control signal is a switch signal (dry contact), an open contact is considered high level, and a closed contact is considered low level. When the signal is high level, the pump will run clockwise. When the signal is low level, the pump will run counterclockwise.

Low Level Signal for CW: When the control signal is a logic level signal, a voltage of 5V to 24V is considered high level. When the control signal is a switch signal (dry contact), an open contact is considered high level, and a closed contact is considered low level. When the signal is high level, the pump will run counterclockwise. When the signal is low level, the pump will run clockwise.

Setting steps:

1. When the pump is on standby and displaying any run screen, press and hold the Function Key



enter the Control Mode setting screen. Then tap Enter Key six times to enter the Running Direction Trigger Mode setting Screen.

2. Tap Increase and Decrease Keys

to set the trigger mode.

- to confirm and save the setting, and display the next parameter setting screen 3. Tap the Enter Kev (Pump Address).
- to return to the previous 4. If no other parameters need to be set, press and hold the Function Key run screen. The pump only can start running when displaying run screen.

3.2.3.8 Pump Address for the Communication Control

When using communication commands to set the parameters and control the pump, the pump address can be set to 1-30.





Pump Address 30

Setting steps:



1. When the pump is on standby and displaying any run screen, press and hold the Function Key



seven times to enter the Pump Address enter the Control Mode setting screen. Then tap Enter Key setting screen.

2. Tap Increase and Decrease Keys

to set the pump address.

- to confirm and save the setting, and display the next parameter setting screen 3. Tap the Enter Key (Baud Rate).
- to return to the previous 4. If no other parameters need to be set, press and hold the Function Key run screen. The pump only can start running when displaying run screen.

3.2.3.9 Baud Rate for the Communication Control

When using communication commands to set the parameters and control the pump, the baud rate can be set to 1200bps (displayed as 1.2), 2400bps (displayed as 2.4), 4800bps (displayed as 4.8), 9600bps (displayed as 9.6), 19200 (19.2), 38400bps (displayed as 384), 115200bps (displayed as 115).



Baud rate is 1200bps



Baud rate is 115200bps

Setting steps:

1. When the pump is on standby and displaying any run screen, press and hold the Function Key

enter the Control Mode setting screen. Then tap Enter Key eight times to enter the Baud Rate setting screen.

- 2. Tap Increase and Decrease Keys to set the baud rate.
- 3. Tap the Enter Key to confirm and save the setting, and display the next parameter setting screen (Parity).
- 4. If no other parameters need to be set, press and hold the Function Key \checkmark to return to the previous run screen. The pump only can start running when displaying run screen.

3.2.3.10 Parity for the Communication Control

When using communication commands to set the parameters and control the pump, the parity can be set to NON, ODD, EVEN.





1. When the pump is on standby and displaying any run screen, press and hold the Function Key



enter the Control Mode setting screen. Then tap Enter Key ¹ nine times to enter the Parity setting

Setting steps:

screen. Or tap Function Key

once to enter the Parity setting screen.

- 2. Tap Increase and Decrease Keys to set the parity.
- 3. Tap the Enter Key to confirm and save the setting, and display the next parameter setting screen (Control Mode).
- 4. If no other parameters need to be set, press and hold the Function Key 🖤 to return to the previous

run screen. The pump only can start running when displaying run screen.

3.3 Flow Rate Calibration

Gx00-1L series pump can display flow rate and pumped volume, and operate according to the preset target volume. To ensure the displayed flow rate and volume match the actual values, and the pumped volume meets the set target volume, it is necessary to calibrate the flow rate before the pump starts official operation, or after replacing the pump head or tubing.

The flow rate can be calibrated by using the actual pumped volume or actual flow rate. And the calibration can be performed in keypad control mode or communication control mode. External control mode does not support the calibration function.

3.3.1 Flow Rate Calibration in Keypad Control Mode

3.3.1.1 Calibration with Volume



Calibration Screen with Volume

Calibration steps:

1. When the pump displays any run screen, tap the Function Key

to display the calibration screen.

- 2. Press and Hold the Unit Switch Key to change the unit to a volume unit.
- 3. Start the pump, and the pump will begin timing automatically.
- 4. Stop the pump at any time. (If the pump was set with target time in Time Screen or target volume in Volume Screen, the pump will stop automatically when the target time is reached or the target volume is reached.)

Tap the Increase and Decrease Keys 5.

pumped volume during the running time.



to calibrate, and the LCD will show "CAL" for 2 seconds.

Note:

6.

- 1. After calibration, the pump speed will adjust to match the flow rate set before calibration.
- 2. After calibration, if the pump speed needs to be set higher than the max speed (the max speed of G100-1L is 150rpm, the max speed of G300-1L is 350rpm, the max speed of G600-1L is 650rpm.), the LCD will show "err" for 2 seconds, and the pump speed will adjust to the max speed, and the flow rate will change to the max flow rate at the max speed. This indicates that the set flow rate exceeds the capacity of the loaded tubing.

3.3.1.2 Calibration with Flow Rate

Tap the Enter Key



Calibration Screen with Flow Rate

If the pump is working in timed transfer mode or volume transfer mode, calibrate the flow rate as below steps:

to display the calibration screen. When the pump displays any run screen, tap the Function Key 1.

լուն to change the unit to a flow rate unit. Press and Hold the Unit Switch Kev 2.

- 3. Start the pump, and the pump will begin timing automatically.
- Stop the pump at any time, or wait for the pump to stop automatically when the target time is reached or 4. the target volume is reached.
- to change the three-digit value according to the actual Tap the Increase and Decrease Keys 5. flow rate corresponding to the current pump speed.



to calibrate, and the LCD will show "CAL" for about 2 seconds. Tap the Enter Key

Note:

6.

- 1. After calibration, the pump speed will adjust to match the flow rate set before calibration.
- 2. After calibration, if the pump speed needs to be set higher than the max speed (the max speed of G100-1L is 150rpm, the max speed of G300-1L is 350rpm, the max speed of G600-1L is 650rpm.), the LCD will show "err" for 2 seconds, and the pump speed will adjust to the max speed, and the flow rate will change to the max flow rate at the max speed. This indicates that the set flow rate exceeds the capacity of the loaded tubing.

If the pump is working in continuous transfer mode, calibrate the flow rate as below steps.

- **\$** to display the calibration screen. When the pump displays any run screen, tap the Function Key 1.
- to change the unit to a flow rate unit. Press and Hold the Unit Switch Key 2.
- 3. Start the pump or not.
- to change the three-digit value according to the actual Tap the Increase and Decrease Keys \downarrow 4. flow rate corresponding to the current pump speed.



to calibrate, and the LCD will show "CAL" for about 2 seconds. Tap the Enter Kev

Note:

5.

- 1. After calibration, the pump speed will adjust to match the flow rate set before calibration.
- 2. After calibration, if the pump speed needs to be set higher than the max speed (the max speed of G100-1L is 150rpm, the max speed of G300-1L is 350rpm, the max speed of G600-1L is 650rpm.), the LCD will show "err" for 2 seconds, and the pump speed will adjust to the max speed, and the flow rate will change to the max flow rate at the max speed. This indicates that the set flow rate exceeds the capacity of the loaded tubing.

3.3.2 Flow Rate Calibration in Communication Control Mode

When the pump is working in Communication Control Mode, the flow rate can be calibrated by receiving the actual pumped volume during the running time, or the actual flow rate. Please refer to Appendix 2 for the Modbus RTU protocol for the register addresses of calibration parameter, calibration value, and calibration unit.

Note:

- 1. After calibration, the pump speed will be changed to the value corresponding to the set flow rate. And the flow rate will remain the same as the one set before calibration.
- 2. After calibration, if the pump speed needs to be changed to a value higher than the max speed (the max speed of G100-1L is 150rpm, the max speed of G300-1L is 350rpm, the max speed of G600-1L is 650rpm.), the LCD will show "err" for 2 seconds, and the pump speed will be changed to the max speed, and the flow rate will be changed to the value corresponding to the max speed. This means the set flow rate is out of the loaded tubing capacity.

4 Definition of the External Control Interface



PIN	Functional Definition	Description	Wiring Instruction			
1	START/STOP	Start/stop signal input	(1) Using a voltage signal or frequency signal to control the pump speed:			
2	CW/CCW	Running direction signal input, CW: clockwise, CCW: counterclockwise				
3	V/F	Pump speed control signal input: Voltage signal: 0-5V or 0-10V Frequency signal: 0-10kHz	(2) Using a current signal to control the pump speed:			
4	mA	Pump speed control signal input: 4-20mA				
5	GND	Ground	START CW/ 4-20mA /STOP CCW 5/8			
6	А	RS485+	 (3) Using an external signal (e.g. foot pedals) to control the start/stop of the pump in keypad control mode: 			
7	В	RS485-	START /STOP 5/8			
8	GND	Ground	 (4) To control the start/stop, and running direction, Pin 1 and Pin 2 can be connected to the 5 to 24V voltage signal or dry contact switch signal 			
9	PE	Protecting Earth	Start/Stop Control: CW/CCW Control: 5/8 1 5/8 1 5/8 1 5/8 2 5/8 2 5V-24V Dry Contact 5V-24V Dry Contact 5V-24V Dry Contact			

5 Appendix

Appendix 1: Parameter Default Setting

Set all parameters to factory default setting: During the power-up process, press and hold the function key

with the LCD displays "---". Then tap the "+ or -" button to display "FFF", and tap the confirmation key

to reset all the parameters to their default settings.

Parameter	Factory Default Setting	Remark		
Pump speed on the Pump Speed	G100-1L: 100rpm	Pump will display the Pump Speed Run Screen when		
Run Screen	G300-1L: 300rpm	nowered up for the first time		
	G600-1L: 600rpm			
	G100-1L: 500 mL/min			
Flow rate on the Flow Rate Screen	G300-1L: 1.5L/min			
	G600-1L: 3L/min			
Target time on the Time Screen	0: 00	If the target time and target volume are 0, the work mode		
Targe volume on the Volume Screen	0.00uL	is continuous transfer mode.		
Pumped volume on Calibration	0.00.01	The conversion factor between flow rate and pump speed		
Screen	0.000	is the factory default setting.		
Control mode	INT	INT means keypad control mode		
Powered up status	Stop			
Keypad lock	Never lock			
External pump speed control signal	0-5V			
	G100-1L: 150 rpm			
Max speed for external control	G300-1L: 350 rpm			
	G600-1L: 650 rpm			
	The default signal range			
	corresponds to pump			
The relationship between en	speed range: 0rpm to the	The valid signal range and corresponding appending		
automatic signal and pump around	max speed:	The value signal range and corresponding speed range can		
external signal and pump speed	G100-1L: 150 rpm	be set through communication commands.		
	G300-1L: 350 rpm			
	G600-1L: 650 rpm			
Start/stop control trigger mode	Falling edge trigger			
Running direction control trigger	Falling edge trigger			
mode				
Pump address	01			
Baud rate	115200bps			
Parity	NON			
Stop bit	1	The communication stop bit cannot be set		

Refer to Appendix 2 for the factory default setting for other parameters in the Modbus RTU protocol.

Appendix 2: Modbus RTU Protocol

To control the pump with communication commands, set the control mode to COM through the keypad or by setting the register address of 0x0020 to 1. The pump address, baud rate, parity and stop bit should be the same as the pump setting.

The default settings are: pump address 1, baud rate 115200bps, no parity, 1 stop bit.

Note: some external signal control parameters can also be set through communication commands. The control mode should be COM before setting those external control parameters.

Parameter	Register	Data	Read/	Be saved after	Factory setting	Description
Vallable	4441035	type	WIIte	power off?		
Stort/stop	0×0001	wint 16		N	0	0: Stop,
Start/stop	0x0001	um_ro	r./ vv	IN	0	1: Start
Run at the full		_				0: Run at the normal speed
speed	0x0006	uint_16	R/W	N	0	1: Run at the full speed
Pump address	0x0010	uint_16	R/W	Y	1	Setting range: 1-30
						0: 1200bps
						1: 2400bps
						2: 4800bps
Baud rate	0x0011	uint_16	R/W	Y	6	3: 9600bps
						4: 19200bps
						5: 38400 bps
						6: 115200bps
						0: Non,
Parity	0x0012	uint_16	R/W	Y	0	1: Odd,
						2: Even
						1: enable the communication control
						mode
						To enable the communication control
Enable						mode, set this register to 1 while the
communication	0x0020	uint 16	R/W	Y	0	pump is in keypad control mode or
control mode		unit_ro			U U	external signal control mode.
control mode						Note that keypad control mode and
						external signal control mode can only
						be set via the keypad, not through
						communication commands.
						0: do not lock the keypad,
Keypad lock	0x0030	uint 16	R/W	Y	0	5 - 60: 5s-60s, the keypad will lock
					-	when it is not operated for the set
						time

External signal control setting: Start/stop trigger mode	0x0031	uint_16	R/W	Y	2	 bit1: bit0 00: Level trigger. Open contact or high level, pump runs. Closed contact or low level, pump stops. 01: Level trigger. Closed contact or low level, pump runs. Open contact or high level, pump stops. 10: Pulse trigger. Falling edge signal switch the status between start and stop 11: Pulse trigger. Rising edge signal switch the status between start and top
External signal control setting: Running direction trigger mode	0x0032	uint_16	R/W	Y	2	 bit1: bit0 00: Level trigger. Open contact or high level, pump runs CW. Closed contact or low level, pump runs CCW. 01: Level trigger. Closed contact or low level, pump runs CW. Open contact or high level, pump runs CCW. 10: Pulse trigger. Falling edge signal switch the direction between CW and CCW 11: Pulse trigger. Rising edge signal switch the direction between CW and CCW.
External signal control setting: Pump speed control signal type	0x0033	uint_16	R/W	Y	0	0: 0-5V, 1: 0-10V, 2: 4-20mA, 3: 0-10KHz
External signal control setting: Max speed for external signal control	0x0034	uint_16	R/W	Y	G100-1L: 150 rpm G300-1L: 350 rpm G600-1L: 650 rpm	(min speed set in 0x0035 + 1rpm) to max speed Max speed: G100-1L 150rpm, G300- 1L 350rpm, G600-1L 650rpm Data unit: 0.01rpm
External signal control setting: Min speed for external signal control	0x0035	uint_16	R/W	Y	0	0rpm to (max speed set in 0x0034 - 1rpm) Data unit: 0.01rpm
External signal control setting: Min input for the 0-5V speed control	0x0036	uint_16	R/W	Y	0	0V to (max input set in 0x0037 -1V) Data unit: 0.01V

External signal control setting: Max input for the 0-5V speed control	0x0037	uint_16	R/W	Y	5V	(min input set in 0x0036 +1V) to 5V Data unit: 0.01V
External signal control setting: Min input for the 0-10V speed control	0x0038	uint_16	R/W	Y	0	0V to (max input set in 0x0039 -1V) Data unit: 0.01V
External signal control setting: Max input for the 0-10V speed control	0x0039	uint_16	R/W	Y	10V	(min input set in 0x0038 +1V) to 10V Data unit: 0.01V
External signal control setting: Min input for the 4-20mA speed control	0x003A	uint_16	R/W	Y	4mA	4mA to (max input set in 0x003B - 1.6mA) Data unit: 0.01mA
External signal control setting: Max input for the 4-20mA speed control	0x003B	uint_16	R/W	Y	20mA	(min input set in 0x003A +1.6mA) to 20mA Data unit: 0.01mA
External signal control setting: Min input for the 0-10kHz speed control	0x003C	uint_16	R/W	Y	0	0 to (max input set in 0x003D-1KHz) Data unit: 1Hz
External signal control setting: Max input for the 0-10kHz speed control	0x003D	uint_16	R/W	Y	10KHz	(min input set in 0x003C +1kHz) to 10kHz Data unit: 1Hz
Running direction	0x0060	uint_16	R/W	Y	1	0: counterclockwise, 1: clockwise
Target volume in volume transfer mode	0x0063	uint_16	R/W	Y	0	0 to 9999, max 999L
Unit for the target volume	0x0064	uint_16	R/W	Y	101	99: 100nL, 100: 1uL, 101: 10uL, 102: 100uL, 103: 1mL, 104: 10mL, 105: 100mL,

						106: 1L
Target time in						
timed transfer	0x0065	uint_16	R/W	Y	0	0 to 9999, max 999 minutes
mode						
						100: 1s,
Unit for the target	0x0066	uint 16	R/W	Y	100	101: 0.1min,
time						102: 1min
						0 to 9999
					G100-1L: 500	G100 flow rate range: 0-750mL/min
Flow rate	0x0067	uint_16	R/W	Y	G300-1L: 150	G300 flow rate range: 0-1 751 /min
					G600-1L: 300	G600 flow rate range: 0-3 251 /min
						99: 100nL/min
						100: 1ul /min
						101: 10ul /min
Linit for the flow					G100-1L: 103	102: 100ul /min
rate	0x0068	uint_16	R/W	Y	G300-1L: 104	102: 1ml /min
late					G600-1L: 104	104: 10ml (min
					G100-1L: 100	0 to 9999
Pump speed	0x0069	uint_16	R/W	Y	G300-1L: 300	G100 speed range: 0-150rpm
					G600-1L: 600	G300 speed range: 0-350rpm
						G600 speed range: 0-650 rpm
Unit for the pump			-			98: 0.01rpm
speed	0x006A	uint_16	R/W	Y	100	99: 0.1rpm
						100: 1rpm
						0: calibrate the pump with pumped
Calibration						volume
variable	0x00BF	uint_16	R/W	N	2	2: calibrate the pump with actual flow
						rate corresponding to the current
						speed
Actual value for					G100-1L: 500	
calibration	0x00C0	uint_16	R/W	N	G300-1L: 150	0 to 9999
					G600-1L: 300	
						Pumped volume:
						99: 100nL,
					G100-11 103	100: 1uL,
Unit for the actual	0x00C1	uint 16	R/W	N	G300-11 · 104	101: 10uL,
value					G600-11 · 104	102: 100uL,
					0000-12.104	103: 1mL,
						104: 10mL,
						105: 100mL,

						106: 1L
						Flow rate:
						97: 1nL/min,
						98: 10nL/min,
						99: 100nL/min,
						100: 1uL/min,
						101: 10uL/min,
						102: 100uL/min,
						103: 1mL/min,
						104: 10mL/min,
						105: 100mL/min,
						106: 1L/min
						Error Code:
						BIT0: abnormal status of external
						signal control
						BIT1: abnormal status of data reading
						BIT2: abnormal status of data saving
						BIT3: the target time or volume is too
System Status	0x0100	uint 16	R	N	0	small to run the pump
(High Letter)	0,0100	unit_10			U U	BIT8: the running time exceeds
						999minutes according to the target
						volume and flow rate
						BIT9: the volume exceeds 999L
						according to the target time and flow
						rate
						BIT10: Calibration out of range

Appendix 3: Longer Pump OEM Communication Protocol

To control the pump with communication commands, set the control mode to COM through the keypad. The pump address, baud rate, parity and stop bit need to be the same as the pump setting.

The default settings are: pump address 1, baud rate 115200bps, no parity, 1 stop bit.

- 1. Frame format (Character format): 1start + 8data + 1parity + 1stop
 - 1 start bit
 - 8 data bits
 - 1 parity bit (non, odd or even)
 - 1 stop bit
- 2. Message format: flag+addr+len+pdu+fcs.
 - flag: E9H, the message head. When sending the message, the data E8H after message head will be replaced with E8H 00H, and E9H after message head will be replaced with E8H 01H. When receiving the message, the data E8H 00H after message head will be reverted to E8H, and E8H 01H after message head will be reverted to E9H. (Note: if E8 00 replaced E8 or E8 01 replaced E9, E8 00 or E8 01 will be regarded as one byte, no influence on the length of pdu.)
 addr: one byte, pump address, 1-30. 31 is broadcast address.
 - len: one byte, length of **pdu**.

fcs: one byte, XOR of addr, len , pdu.

3. **pdu** format: application layer code format

3.1 Set the pump running parameters

Send to the pump:



Pump responds:



- WJ: 2 bytes, using ASCII code, to indicate that this command is used to set pump running parameters. ASCII code of W is 57H, ASCII code of J is 4AH.
- Speed: 2 bytes. hexadecimal number, most significant byte first. The data unit is 0.01rpm. Example: 2710H means 100rpm.
- Full speed and start/stop:

Bit0: 1 means pump runs, 0 means pump stops.

Bit1: 1 means full speed, 0 means normal speed.

• Direction:

Bit0: 1 means CW, 0 means CCW.

Note: When setting the running parameter, the **addr** in message can be pump address (1-30) or broadcast address 31. All pumps will operate according to the same command without response when using broadcast address.

3.2 Read the pump running status

Send to the pump:

RJ

Pump responds:



Speed (2 bytes)

Full speed and start/stop (1 byte)

Direction (1 byte)

- RJ: 2 bytes, using ASCII code, to indicate that this call is used to read pump running status. ASCII code of R is 52H, ASCII code of J is 4AH.
- Refer to para. 3.1 for instructions of speed, full speed, start/stop, direction

Note: When reading the running state, the addr in the message only can be pump address (1-30).

3.3 Set the flow rate parameters

Send to the pump:



Pump responds:

WL

- WL: 2 bytes, using ASCII code, to indicate that this command is used to set the flow rate parameters. ASCII code of W is 57H, ASCII code of L is 4CH.
- Flow rate: 4 bytes. hexadecimal number, most significant byte first. The data unit is 1nL/min. 1L=10³ mL=10⁶ µL=10⁹ nL Example: 5F5 E100H means 100mL/min.
- Refer to para. 3.1 for instructions of full speed, start/stop, direction

Note: When setting the flow rate parameters, the **addr** in message can be pump address (1-30) or broadcast address 31. All pumps will operate according to the same command without response when using broadcast address.

3.4 Read the flow rate parameters

Send to the pump:



Pump responds:



Full speed and start/stop (1 byte)

Direction (1 byte)

- RL: 2 bytes, using ASCII code, to indicate that this command is used to read flow rate parameters. ASCII code of R is 52H, ASCII code of L is 4CH.
- Flow rate: 4 bytes. hexadecimal number, most significant byte first. The data unit is 1nL/min. 1L=10³ mL=10⁶ µL=10⁹ nL Example: 5F5 E100H means 100mL/min.
- Refer to para. 3.1 for instructions of full speed, start/stop, direction

Note: When reading the running state, the addr in the message only can be pump address (1-30).

3.5 Set the target time for timed transfer mode

Send to the pump:



Pump responds:

WM

- WM: 2 bytes, using ASCII code, to indicate that this command is used to set the target time. ASCII code of W is 57H, ASCII code of M is 4DH.
- Unit of the time: 100 means 1s, 101 means 10s, 102 means 1min.
- Target time: 2 bytes, hexadecimal number, most significant byte first.
 When the unit of the time is 1s, the target time can be set to 0-600, means 0-600s
 When the unit of the time is 10s, the target time can be set to 0-600, means 0-6000s
 When the unit of the time is 1min, the target time can be set to 0-999, means 0-999mins
- Refer to para. 3.1 for instructions of full speed, start/stop, direction

Note: When setting the target time parameter, the **addr** in message can be pump address (1-30) or broadcast address 31. All pumps will operate according to the same command without response when using broadcast

3.6 Read the target time parameters

Send to the pump:



Pump responds:

RM	Target time (2 bytes)	The unit of the time (1 byte)	Full speed and start/stop (1 byte)	Direction (1 byte)
----	--------------------------	-------------------------------	------------------------------------	--------------------

- RM: 2 bytes, using ASCII code, to indicate that this command is used to read the target time parameters. ASCII code of R is 52H, ASCII code of M is 4DH.
- Unit of the time: 100 means 1s, 101 means 10s, 102 means 1min.
- Target time: 2 bytes, hexadecimal number, most significant byte first.
 When the unit of the time is 1s, the target time can be set to 0-600, means 0-600s
 When the unit of the time is 10s, the target time can be set to 0-600, means 0-6000s
 When the unit of the time is 1min, the target time can be set to 0-999, means 0-999mins
- Refer to para. 3.1 for instructions of full speed, start/stop, direction

Note: When reading the target time, the **addr** in the message only can be pump address (1-30).

3.7 Set the target volume for volume transfer mode

Send to the pump:

WV	Target volume (2 bytes)	The unit of the volume (1 byte)	Full speed and start/stop (1 byte)	Direction (1 byte)
----	----------------------------	---------------------------------	---------------------------------------	--------------------

Pump responds:

WV

- WV: 2 bytes, using ASCII code, to indicate that this command is used to set the target volume. ASCII code of W is 57H, ASCII code of V is 56H.
- Target volume: 2 bytes, hexadecimal number, most significant byte first. Setting range: 0-999
- Unit of the volume: 1byte, setting range 98-106

```
98-0.01 uL, \ 99-0.1 uL, \ 100-1 uL, \ 101-0.01 mL, \ 102-0.1 mL, \ 103-1 mL, \ 104-0.01 L, \ 105-0.1 L, \
```

106 - 1L

• Refer to para. 3.1 for instructions of full speed, start/stop, direction

Note: When setting the target volume parameter, the **addr** in message can be pump address (1-30) or broadcast address 31. All pumps will operate according to the same command without response when using broadcast address.

3.8 Read the target volume parameter

Send to the pump:

RV

Pump responds:

RV	Target volume (2 bytes)	The unit of the volume (1 byte)	Full speed and start/stop (1 byte)	Direction (1 byte)
----	----------------------------	---------------------------------	------------------------------------	--------------------

- RV: 2 bytes, using ASCII code, to indicate that this command is used to read the target volume. ASCII code of R is 52H, ASCII code of V is 56H.
- Target volume: 2 bytes, hexadecimal number, most significant byte first. Setting range: 0-999
- Unit of the volume: 1byte, setting range 98-106

 $98-0.01 uL, \ 99-0.1 uL, \ 100-1 uL, \ 101-0.01 mL, \ 102-0.1 mL, \ 103-1 mL, \ 104-0.01 L, \ 105-0.1 L, \$

106 - 1L

• Refer to para. 3.1 for instructions of full speed, start/stop, direction

Note: When reading the target volume, the **addr** in the message only can be pump address (1-30).

3.9 Set communication parameters:

Send to the pump:



Pump responds:

WID

- WID: 3 bytes, using ASCII code, to indicate that this command is used to set communication parameters. ASCII code of W is 57H, ASCII code of I is 49H, ASCII code of D is 44H.
- Pump address: 1 bytes, hexadecimal number, setting range: 0-30
- Baud rate: 2 bytes, hexadecimal number, most significant byte first. Setting range:1-7

1-1200bps, 2-2400bps, 3-4800bps, 4-9600bps,

5-19200bps, 6-38400bps, 7-115200bps

- Parity: 1 byte, 1-Non, 2-Odd, 3-Even
- Stop bit: 1 byte, only can be set to 1.

Note: When setting the communication parameter, the **addr** in the message can be pump address (1-30) or broadcast address 31. All pumps will operate according to the same command without response when using broadcast address.

3.10 Calibrate the pump with actual flow rate

Send to the pump:

CL

The actual flow rate (4 bytes)

Pump responds:



- CL: 2 bytes, using ASCII code, to indicate that this command is used to calibrate the pump with actual flow rate corresponding to the current speed. ASCII code of C is 43H, ASCII code of L is 4CH.
- Actual flow rate: 4 bytes. hexadecimal number, most significant byte first. The data unit is 1nL/min.
 1L=10³ mL=10⁶ µL=10⁹ nL Example: 5F5 E100H means 100mL/min.
- 3.11 Calibrate the pump with pumped volume:

Send to the pump:



Pump responds:

CV

- CV: 2 bytes, using ASCII code, to indicate that this command is used to calibrate the pump with actual pumped volume. ASCII code of C is 43H, ASCII code of V is 56H.
- Actual pumped volume: 2 bytes. hexadecimal number, most significant byte first. Setting range: 0-999.
- Unit for the volume: 1 byte, setting range: 98-106

98 - 0.01uL, 99 - 0.1uL, 100 - 1uL, 101 - 0.01mL, 102 - 0.1mL,

103 - 1mL, 104 - 0.01L, 105 - 0.1L, 106 - 1L

4 Examples:

(1) To set G100-1L (addr: 01) to run clockwise, the speed is 100rpm. The command string is as follows E9 01 06 57 4A 27 10 01 01 2E

(2) To set G100-1L (addr: 01) to run clockwise, and the flow rate is 100mL/min. The command string is as follows

E9 01 08 57 4C 05 F5 E1 00 01 01 03

Appendix 4: Alarm Display

Error Codes	Description	Remark
E01	The communication with the external	Contact after-sales service
	control board is abnormal	
Err	Calibration out of range	Check the calibration value input:
		1. If the pump speed after calibration will exceed the max speed, the
		displayed speed will adjust to the max speed, and the flow rate will
		change to the max flow rate corresponding to the max speed.
		2. If the flow rate after calibration will exceed the pump capacity, the
		calibration will not be performed.
CAL	The calibration was performed	
	successfully	