

OXY-1 ST (trace) & OXY-4 ST (trace) (G3)

TRANSMITTERS

[Instruction Manual](#)



OXY-1 ST (trace) & OXY-4 ST (trace) (G3)

Specification:

Fiber optic (trace) oxygen meter for use with oxygen sensor spots, flow-through cells, dipping probes & microsensors

Software:

PreSens Measurement Studio (PMS2)

Document filename: IM_OXY-1 -4 ST_dv3

All rights reserved. No parts of this work may be reproduced in any form or by any means - graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems - without the written permission of the publisher.

Products that are referred to in this document may be either trademarks and/or registered trademarks of the respective owners. The publisher and the author make no claim to these trademarks.

While every precaution has been taken in the preparation of this document, the publisher and the author assume no responsibility for errors or omissions, or for damages resulting from the use of information contained in this document or from the use of programs and source code that may accompany it. In no event shall the publisher and the author be liable for any loss of profit or any other commercial damage caused or alleged to have been caused directly or indirectly by this document.

Specifications may change without prior notice.

Manufacturer

PreSens

Precision Sensing GmbH

Am BioPark 11

93053 Regensburg, Germany

Phone +49 941 94272100

Fax +49 941 94272111

info@PreSens.de

www.PreSens.de

Table of Contents

1	Preface.....	7
2	Safety Notes.....	8
3	Description of the OXY-1/-4 ST (trace).....	9
4	Set-Up	10
4.1	Connecting the Devices	10
4.2	Important Considerations for USB Handling	12
5	Technical Data	13
6	Operational Notes	14
6.1	Optical Output	14
6.2	Maintenance.....	14
6.3	Service.....	14
7	Concluding Remarks	15

1 Preface

You have chosen a new, innovative technology for measuring oxygen.

The OXY-1/-4 ST (trace) is a compact fiber optic oxygen meter. It is based on a novel technology, which creates very stable, internally referenced measured values. This allows a more flexible use of oxygen sensors in various fields of interest.

Optical oxygen sensors (also called optrodes) have several important features:

- They are small.
- Their signal does not depend on the flow rate of the sample.
- They can be physically divided from the measuring system which allows a non-invasive measurement.
- They can be used in disposables.
- Microsensors allow measurements with high spatial resolution whenever this is required.

Therefore, they are ideally suited for the examination of small sample volumes, for highly parallelized measurements in disposables, and for biotechnological applications. A set of different oxygen mini- and microsensors, flow-through cells and non-invasive sensors is available to make sure you have the sensor which matches your application.

Please feel free to contact our service team to find the best solution for your application.

Your PreSens Team

PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE WORKING WITH THIS DEVICE. WHEN DISREGARDING THESE INSTRUCTIONS THE SAFETY OF THE DEVICE CAN BE IMPAIRED.

2 Safety Notes

! It is the customer's responsibility to validate the sensor and transmitter under end-user conditions according to safety precautions of the application to ensure that the use of the sensor is safe and suitable for the intended purpose.

PreSens is explicitly not liable for direct or indirect losses caused by the application of these measurement systems. In particular it has to be considered that malfunctions can occur due to the naturally limited lifetime of the sensor depending on the respective application. The set-up of backup measurement stations is recommended when using the sensors in critical applications to avoid consequential losses. It is the customer's responsibility to install a suitable safety system in the event of sensor failure.

3 Description of the OXY-1/-4 ST (trace)



Fig. 1 OXY-4 ST (G3) and OXY-1 ST

These oxygen meters can be used with microsensors (230 μm fiber) in different designs, non-invasive sensors, dipping probes and flow-through cells (1 mm fiber). OXY-1/-4 ST are compatible with sensor type PSt7 (detection limit 15 ppb, 0 - 100 % oxygen), and the trace versions are additionally compatible with sensor type PSt8 (detection limit 3 ppb, 0 - 10 % oxygen). OXY-1/-4 ST devices are USB powered and do not need an extra power adapter. All OXY-1/-4 ST have temperature compensation for each individual channel, so most precise measurements in environments with changing temperature can be taken. It is also possible to group several channels, and use the measurement values of one temperature sensor for compensation on all grouped sensors. The oxygen meters are operated with the PreSens Measurement Studio 2 software, which enables simultaneous control of several devices. With numerous features and additional pressure and salinity compensation, the software makes the OXY-1/-4 ST suitable for numerous applications.

! Several single- and multichannel devices can be connected to the PreSens Measurement Studio 2 software, so you can measure with as many sensors simultaneously as required.

4 Set-Up

4.1 Connecting the Devices

1. Remove the protective cap from the male plug on the optical fiber and insert it in the ST connection of the OXY-1/-4 ST. The ST plug has to be inserted and slightly turned clockwise to fasten it. Be careful not to snap off the optical fiber; best hold the fiber between forefinger and thumb at the bayonet lock of the male plug.



Fig. 1 Connecting the ST plug to the OXY-1 ST or OXY-4 ST

For further information on sensor handling and calibration, please refer to the respective sensor instruction manual (www.presens.de/manuals).

2. Connect the USB cable to the connection on the OXY-1 ST and to an USB port of your PC / notebook.



Fig. 2 Connecting the USB cable

3. If required, you can now connect a Pt100 temperature sensor to the device.



Fig. 3 Connecting the temperature sensor to OXY-1 ST

There is a red mark on the temperature sensor connection of the OXY-1 ST as well as on the temperature sensor plug. Match those two marks before inserting the temperature sensor plug into the connection on the device, else the plug might get damaged.



Fig. 4 Connecting the temperature sensor to OXY-4 ST

For connecting the temperature sensor to an OXY-4 ST the red mark on the temperature sensor plug has to point downward. Turn the plug, so the red mark points down and carefully insert it into the connector on the OXY-4 ST.

4. Now start the PreSens Measurement Studio (PMS2) software on your PC / notebook. (The OXY-1 ST always has to be connected via USB before the software is started.) After successful initialization the software main screen is displayed and the OXY-1 ST will show in the **Device** section.

For further details on how to use the PreSens Measurement Studio please refer to the software instruction manual.

4.2 Important Considerations for USB Handling

- ! PreSens recommends the use of a dedicated USB 2.0 PCI Card to connect and handle USB PreSens Devices with a Desktop PC.
- ! In order to enhance the system stability avoid the use of USB hubs and connect PreSens Devices directly to your PC / notebook USB Ports.
- ! If possible, disconnect all other USB devices that are not in use, as they may reduce or disturb the USB resources of your PC / notebook.
- ! Docking stations may also reduce or disturb the USB resources of your PC / notebook and therefore affect the correct function of the software.
- ! It is also recommended to disable the Power Saving Settings of your USB Root controller.

5 Technical Data

OPTICAL SENSOR

Oxygen sensor	OXY-1/-4 ST: PSt7 OXY-1/-4 ST trace: PSt7, PSt8
Optical connector	ST
Channels	1 / 4; up to 10 OXY-1/-4 ST (trace) can be controlled simultaneously via PMS2 software
LED peak wavelength	633 nm

TEMPERATURE SENSOR

Potentiometric temperature sensor (Pt100)	Range	0 – 50 °C
	Resolution	± 0.1 °C

POWER SUPPLY

Supply voltage	5 VDC (USB-2.0-Mini-B)
Current / Power	400 mA

DIGITAL INTERFACE

USB interface cable to PC	Cable included
---------------------------	----------------

ENVIRONMENTAL CONDITIONS

Operating temperature	0 °C to 50 °C
Storage temperature	- 20 °C to 70 °C
Relative humidity	0 % to 80 % (non-condensing)

DIMENSIONS / WEIGHT

W x L x H	OXY-1 ST (trace): ca. 35 x 99 (with connectors) x 30 mm
	OXY-4 ST (trace): ca. 50 mm x 164 mm x 85 mm (with connectors)
Weight	OXY-1 ST (trace): 128 g OXY-4 ST (trace): 417 g

6 Operational Notes

6.1 Optical Output

The ST connector is a high precision optical component. Please keep it clean and dry. Always use the rubber cap to close the output when not in use.

6.2 Maintenance

The transmitter is maintenance-free.

The housing should be cleaned with a cloth only. Avoid any moisture entering the housing. Never use benzene, acetone, alcohol or any other organic solvents.

The ST fiber connector of the sensor can be cleaned with lint-free cloth or a cleaning implement for ST connectors only.

6.3 Service

Alignment, rework or repair work may only be carried out by the manufacturer:

**PreSens
Precision Sensing GmbH**

Am BioPark 11
93053 Regensburg
Germany

Phone +49 941 94272100
Fax +49 941 94272111

info@PreSens.de
www.PreSens.de

Please contact our service team in case of any question. We look forward to helping you and are open for any proposition or criticism.

7 Concluding Remarks

Dear Customer,

With this manual, we hope to provide you with an introduction to work with the OXY-1/-4 ST (trace) fiber optic oxygen meter.

This manual does not claim to be complete. We are endeavored to improve and supplement this version.

We are looking forward to your critical review and to any suggestions you may have.

You can find the latest version at www.PreSens.de.

With best regards,

Your PreSens Team



Manufacturer

**PreSens
Precision Sensing GmbH**

Am BioPark 11
93053 Regensburg
Germany

Phone +49 941 94272100
Fax +49 941 94272111

info@PreSens.de
www.PreSens.de