



Hint: Fluidic 1110 can only be used with every second tank interface of Fluidic 219, due to dimensional constraints. Fluidic 234/ 235 are better suited for the use of all tank interfaces in Fluidic 219.

all values without dimensions are in mm

		Tolerance ISO 8015 General Tolerance ISO 2768-mK
		Scale 1:1
Interaction tank – single tank – 4.5 ml w/o cap	Fluidic 1110	Created/changed on: 30.03.2020
	Page (1/1) Size A4	

MF00037_A4