



OVERVIEW

The fluorinated oil from Emulseo comprises a polyfluorinated surfactant (<0.1%) solubilized in Novec™ 7500 oil. It is especially recommended for solubilization of surfactants from Emulseo.

Fluo-Oil 7500 has especially been formulated to ensure stability, biocompatibility, and reproducibility of microfluidic experiments.

BENEFITS

- Biocompatibility
- Reproducibility
- Stability

HOW TO USE

Emulseo's Fluo-Oil 7500 should be used to dilute Emulseo's surfactants neat or already diluted in HFE-7500 oil.

STORAGE

Emulseo's fluorinated oil should be stored at room temperature protected from light in glass vials.

SAFETY NOTES

- Avoid inhalation, ingestion and contact with skin and eyes
- Wear chemical resistant personal protective equipment
- Wear suitable gloves
- Wear suitable protective clothing
- Safety glasses with side-shields

For additional information, see the Material Safety Data Sheet



CONTACT

If you have any queries, please do not hesitate to e-mail us at: contact@emulseo.com