

## PDC-VPE/VPE-2 Basic Vacuum Pump



Pumping speed 50/60 Hz .....	5.1 m <sup>3</sup> /hr
Ultimate total pressure (gas ballast closed) .....	7.5 mTorr (0.01 mbar)
Ultimate total pressure (gas ballast open) .....	150 mTorr (0.2 mbar)
Inlet .....	1/4" Female BSPP
Outlet .....	1/2" Female BSPP
Motor power 50/60 Hz .....	373 W
Noise level (at 3 feet) .....	54 dBA
Weight .....	27 lbs (12 kg)
Dimensions .....	9.75" H x 14" W x 5.5" D
	25 cm H x 36 cm W x 14 cm D