



	ITEM NUMBER: 03033	DRAWN BY: WVDC DRAWN DATE: 01 Jun 2021 PROJECT NO: RD0946.001 EXT. REF.:	TITLE: Chip Serpentine 140x250um, 500mm		
	THIRD ANGLE PROJECTION 	MATERIAL: Glass D263 Bio	DWG NO. 42004744	REVISION 0	SCALE: 2:1
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MICRONIT BV. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MICRONIT BV IS PROHIBITED.	UNLESS OTHERWISE SPECIFIED: - DIMENSIONS ARE IN MILLIMETERS - LINEAR TOLERANCES: ± 0.1 - ANGULAR TOLERANCES: $\pm 5^\circ$	DO NOT SCALE DRAWING	SHEET FORMAT: A4	SHEET 1 OF 1	